

Release Notes 1.3.1

New features

- Printer LP4024, printout in Chinese language available (no settings on printer required).
Note: Printer Firmware LP 4024 Vers:3.0 required

Max 225g Min 1mg d=0.01/0.1mg e=1mg		Administrator - 22.05.2024 16:02	
		系统设置  	
称重参数		模式	协议
数据交换	▲	波特率	9600
串口		奇偶校验	8-无-1
协议	▼	握手协议	无
重量转换		格式化	字符串40
外设设置		打印机	LP4024
通用设置			

Miscellaneous

- Translations revised
- Context help for addon Identification added

Bug fixes

- Add-On Identification: Customizing of soft buttons [ID1..6]
- Add-On Checkweigher: Function of soft button [Capture nominal]